

This belt is manufactured in specially treated Antimicrobial polyester webbing for protection against bacterial infections. It also has a PU Coating to resist the absorption of water.

Recommended washing instructions

To maintain the above properties:

- Hand wash at 40 degrees centigrade using bio or non-bio detergents.
- Do not tumble dry.
- Guaranteed for up to 90 washes.
- Do not use any form of bleach.

Important:

For users where sterile conditions are paramount it is advisable to replace the belt after the recommended number of washes.

This belt is manufactured in specially treated Antimicrobial polyester webbing for protection against bacterial infections. It also has a PU Coating to resist the absorption of water.

Recommended washing instructions

To maintain the above properties:

- Hand wash at 40 degrees centigrade using bio or non-bio detergents.
- Do not tumble dry.
- Guaranteed for up to 90 washes.
- Do not use any form of bleach.

Important:

For users where sterile conditions are paramount it is advisable to replace the belt after the recommended number of washes.