

CRELLING HARNESSSES FOR DISABLED LTD
BELTS AND HARNESSSES FOR SPECIAL NEEDS

WASHING INSTRUCTIONS

- HAND WASH ONLY - In warm water up to 70/80 degrees.
- Use mild detergent/anti-bacterial liquid wash.
- Mild disinfectant can be used.
- DO NOT USE BLEACH
- DO NOT TUMBLE DRY

CRELLING HARNESSSES FOR DISABLED LTD
BELTS AND HARNESSSES FOR SPECIAL NEEDS

WASHING INSTRUCTIONS

- HAND WASH ONLY - In warm water up to 70/80 degrees.
- Use mild detergent/anti-bacterial liquid wash.
- Mild disinfectant can be used.
- DO NOT USE BLEACH
- DO NOT TUMBLE DRY