


Retaining Chair Strap

Some seats may be narrow at the top and the shoulder straps may tend to slip off. Fitting a retaining chair strap horizontally around the seat back - as shown below - will hold the shoulder straps in position.



Cat. Ref: RCS19

The Sketch above shows the Model 19A harness fitted with a RCS19

Crelling Harnesses For Disabled Ltd.

12 Crescent East, Thornton-Cleveleys, Lancashire, FY5 3LJ

Tel: 01253 852298 Fax: 01253 821780 Email: info@crelling.com Website: www.crelling.com